



GENTEC click-style oxygen regulators are available in standard and mini versions, with flow rates of 0-8 LPM, 0-15 LPM, and 0-25 LPM. CGA 540 and CGA 870 inlet connections are available.

These green-anodized aluminum-body regulators have a solid brass core for safety, 1-1/2" diameter UL-listed gauges with gauge protectors, a locking ring on the t-handle to prevent accidental removal, and the FDA and NIOSH-recommended sealing-type washers.

A hose barb fitting is standard on all regulators, while one or two DISS male (with check) power take-offs are available on the standard 0-15 LPM and 0-25 LPM models (in addition to the hose barb).

All GENTEC medical products are backed by our standard 5-year warranty (see warranty statement for details).

284MA-15LY1 Click-Style Regulator



285MA-8LY MINI Click-Style Regulator



## GENTEC Click-style Oxygen Regulator Ordering Information

Model Number	Description
284MA-8LY	Standard CGA 870 regulator, 0-8 LPM, with hose barb
284MA-15LY	Standard CGA 870 regulator, 0-15 LPM, with hose barb
284MA-25LY	Standard CGA 870 regulator, 0-25 LPM, with hose barb
284MA-15LY1	Standard CGA 870 regulator, 0-15 LPM, with hose barb and 1 DISS
284MA-15LY2	Standard CGA 870 regulator, 0-15 LPM, with hose barb and 2 DISS
284MA-25LY1	Standard CGA 870 regulator, 0-25 LPM, with hose barb and 1 DISS
284MA-25LY2	Standard CGA 870 regulator, 0-25 LPM, with hose barb and 2 DISS
285MA-8LY	Mini CGA 870 regulator, 0-8 LPM, with hose barb
285MA-15LY	Mini CGA 870 regulator, 0-15 LPM, with hose barb

NOTE: For CGA 540 inlet connection, drop the "Y" (e.g. 284MA-8L or 285MA-15L).

FLOW RANGE	FLOW SETTINGS (LPM)
0-8 LPM	0 - 0.5 - 1 - 1.5 - 2 - 2.5 - 3 - 4 - 5 - 6 - 7 - 8
0-15 LPM	0 - 0.5 - 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 10 - 15
0-25 LPM	0 - 0.5 - 1 - 2 - 3 - 4 - 5 - 6 - 8 - 10 - 15 - 25

### Specifications

Material: Anodized aluminum body, brass core  
chrome-plated brass t-handle and fittings (hose barb/DISS)

Length: Standard 5-5/8" (14.3 cm)  
Mini 4-5/8" (11.7 cm)

Weight: Standard 12.8 oz. (363 grams)  
Mini 9.6 oz. (272 grams)

Diameter: 1-3/8" (3.5 cm)

Maximum inlet pressure: 3,000 psig (20.7 MPa)

Outlet pressure: 50 psig (345 kPa)

Genstar Technologies Co., Inc.  
4525 Edison Avenue  
Chino, CA 91710 USA  
TEL 909-606-2726  
FAX 909-606-6485  
www.gentechhealthcare.com

