

NEW

DIGITAL SUCTION REGULATORS

Gentec is proud to announce the release of the all-new Digital Suction Regulator. While preserving the look-and-feel and usability of our original 2 1/2" gauge analog version, we provide the added convenience of an easy-to-read digital display (including a digital sweeping needle.) Truly the best of both worlds!

Agion®* Antimicrobial Product Protection

The built-in long-lasting product protection prevents growth of microbes on the surface of the knobs. The antimicrobial product protection is based on non-volatile, naturally occurring silver.



DIGITAL

ANALOG

- Lightweight with durable ABS plastic casing
- Impact-resistant internal components
- Corrosion-free components and can be thoroughly flushed
- On/off cycle settings can be adjusted without opening the case
- Accuracy: +/-1% for digital; +/- 2% for analog
- Intermittent timing cycles: 16 seconds on, 8 seconds off. (timing is adjustable)
- Suction traps and a variety of fittings & adapters are available
- Color-coding available
- 5 year warranty

* Agion's silver antimicrobial has been certified for its material content, recyclability, and manufacturing characteristics. Cradle To Cradle Certified™ products meet established standards for human health, environmental health and recyclability. Cradle To Cradle Certified™ is a certification mark of MBDC.

ORDERING MATRIX

Instructions: Select the suction regulator, adapter, fitting, and color (if needed), then write the model number below.

STEPS:	1 Suction Regulator	2 Adapter	3 Fitting	4 Color (Optional)	Model Number
Example:	881VR-160	- OH	- T	- YL	= 881VR-160-OH-T-YL
Please Write Your Selection:	_____	+ _____	+ _____	+ _____	= _____

1 SUCTION REGULATORS	2 ADAPTERS	3 FITTINGS	4 COLORS (OPTIONAL)	
<p>881VR-160D 881VR-160 Continuous / Intermittent 0-160</p>	<p>1/8" NPT Female 00</p>	<p>1/8" NPT Female A</p>	<ul style="list-style-type: none"> Blank (White) YL (Yellow) BL (Blue) RD (Red) PL (Purple) PC (Peach) PK (Pink) GY (Gray) GN (Green) BB (Baby Blue) 	
<p>881VR-300D 881VR-300 Continuous / Intermittent 0-300</p>	<p>Ohio Adapter 3111U-VAC-NT2 OH</p>	<p>Chemetron Adapter 3103-VAC-NT2 CH</p>		
<p>882VR-300D 882VR-300 Continuous Only 0-300</p>	<p>DISS Hand-Tight Adapter 3222U-VAC-NT2 DH</p>	<p>Trap Fitting 881VR-14K B</p>		
<p>882VR-760D 882VR-760 Continuous Only 0-760</p>	<p>BS: 5682-1998 Adapter 3117-VAC-NT2 BS</p>	<p>DISS Nut-Nipple Adapter 3102-VAC-NT2 DS</p>		
	<p>DIN: 13260-2 Adapter 3116-VAC-NT2 GS</p>	<p>AFNOR: NFS 90-115 Adapter 3105E-VAC-NT2 FS</p>	<p>Complete Suction Trap w/ Trap Fitting 880VT-5 T</p>	
			<p>Complete Suction Trap w/ DISS 880VT-6 TD</p>	<p>Complete Suction Trap w/ DISS Hand-Tight 880VT-7 TH</p>

ACCESSORIES

<p>Suction Trap w/ No Fitting on Top 880VT-0</p>	<p>Suction Trap w/ Male Trap Fitting Only 880VT-1</p>	<p>Suction Trap w/ DISS Female Nut & Nipple 880VT-2</p>	<p>Suction Trap w/ DISS Female Hand-Tight 880VT-3</p>	<p>Suction Trap w/ Tubing Nipple 880VT-4</p>
<p>Vacuum Adapter DISS to Hose Barb DISSM-V-HCHK</p>	<p>Extension Adapter 2-1/2" Long EXTADPT-1/8</p>	<p>Security Cable FMS-2</p>		