



www.gentec.com



EMS GAS FLOW SOLUTIONS

Solutions for Life



EMS GAS FLOW SOLUTIONS

GENTEC® GMAP Series provides a complete gas flow solution for ambulances and emergency care vehicles. Gentec offers gas switchover systems, regulators, fittings, hose assemblies, gas outlets, gas alarms, outlets, flowmeters, suction devices, flowmeters, and EMS accessories. Customized solutions are also available.

CE
0123



1 Medical Regulator

- Regulator for medical gas supply from cylinder
- Chrome-plated brass body
- Sintered inlet filter for longer service life
- CE listed

2 Control Panel

- Continuous monitoring with digital display
- Digital display for ease of monitoring
- Built alarm system, with mute option
- Stainless steel plate, easy to clean

3 Adapters and Fittings

- DISS connection with built-in valve which shuts off gas supply upon disconnection



4 Gas Outlets

- Available in many international standards
- Gas specific outlets
- CE, UL listed

5 Medical Hoses

- Color-coded based on different gases
- Maximum working pressure: 200psi
- Available in custom lengths and configurations

6 Oxygen Therapy

- Dual-sided and elongated scale Flowmeters
- Click-style flowmeters, smallest and lightest in the industry
- Polycarbonate humidifiers bottles, autoclavable

7 Respiratory Therapy

- Oxygen Mask, Nasal Cannula, Breathing Circuits and more emergency equipment products
- CE listed

Ordering Informations

Medical Regulator	
Model Number	Description
195M-CZ	Diaphragm-type Regulator for Ambulance, Chrome-plated Brass Body
197M-CZ	Diaphragm-type Regulator for Ambulance, Matte-chrome Brass Body

Control Panel	
Model Number	Description
GV2-D-DISS-00	Ambulance Control Panel w/ Digital Display and Alarm
CZXT-BOX001	Analog Switchover Control Panel

Adapters and Fittings	
Model Number	Description
T-BFNT2-00	T-Connector
Y-FNT2	Y-Connector
FM-PTO	Power Take-off Connector
FM-IPPB-O	Power Take-off 360-degree Connector
FM197B-20	Rail Clamp
FM197B-20-HB4	Rail Clamp with Hose Barb

Gas Outlets	
Model Number	Description
3851-O	Ohmeda Console Outlet w/ DISS end
3852-O	DISS Console Outlet w/ DISS end
3853-O	Chemetron Console Outlet w/ DISS end
3854-O	Puritan-Bennett Console Outlet w/ DISS end
3900A-11	Trim Plate

* Other international standard adapters are available upon request

Medical Hoses
Available in USA/ISO Standard color, custom length and configurations

Oxygen Therapy	
Model Number	Description
FM197A-15L-XX	Aluminum Body Oxygen Flowmeter, 15LPM
FM197B-15L-XX	Brass Body Oxygen Flowmeter, 15LPM
FM897A-15L-XX	Aluminum Body Oxygen Flowmeter, 15LPM w/ Humidifier
FM897B-15L-XX	Brass Body Oxygen Flowmeter, 15LPM w/ Humidifier
CFM-4L-XX	Click-style Oxygen Flowmeter, 4 LPM
CFM-8L-XX	Click-style Oxygen Flowmeter, 8 LPM
CFM-15L-XX	Click-style Oxygen Flowmeter, 15 LPM
CFM-25L-XX	Click-style Oxygen Flowmeter, 25 LPM
7100R	Polycarbonate Reusable Humidifier, 200ml
7200R	Polycarbonate, Aluminum Top, Reusable Humidifier, 140ml
7300R	Polycarbonate, Reusable Humidifier, 350ml

XX: Inlet adapter (00, OH, DS, DH, CH, PB, FS, GS, BS, JIS, SIS)



Genstar Technologies Company Inc.

4525 Edison Avenue, Chino, CA 91710, USA

Tel: 909-606-2726

Fax: 909-606-6485

www.gentec.com

