



## CLICK-STYLE EMS OXYGEN REGULATORS

GENTEC click-style EMS oxygen regulators are available with either CGA 870 or CGA 540 inlet connections. These durable regulators have a 100% solid brass body, 1-1/2" diameter UL-listed gauge with gauge protector, a locking ring on the t-handle to prevent accidental removal, and the FDA and NIOSH-recommended sealing-type washers. A hose barb fitting is standard on all regulators, while one or two DISS male (with check) fittings are available on all 0-25 LPM models (in addition to the hose barb). GENTEC regulators offer 11 flow positions and large, easy to read flow rates. All GENTEC medical products are backed by our standard 5-year warranty.

### 286MB-25LY2/Left To Right/Barb & 2-DISS



ALL BRASS

### 286MB-25LY/Left to Right/CGA 870/Hose barb only



ALL BRASS

### 286MB-25LY1/Left to Right/CGA 870/Hose barb and 1 DISS



ALL BRASS

# GENTEC Click-style EMS Oxygen Regulator Ordering

Model Number	Description: Yoke Pin Indexed CGA 870 Regulators with Fitting
286MB-25LY2	Standard CGA 870 regulator, 0-25 LPM, with hose barb and 2 DISS
286MB-25LY1	Standard CGA 870 regulator, 0-25 LPM, with hose barb and 1 DISS
286MB-25LY	Standard CGA 870 regulator, 0-25 LPM, with hose barb
Model Number	Description: CGA 540 style regulators with DISS fitting
286MB-25L2	Standard CGA 540 regulator, 0-25 LPM, with hose barb and 2 DISS
286MB-2CK	Standard CGA 540 regulator, 40-60 PSI, with 2 DISS (use on H & M cylinder)

FLOW RANGE	FLOW SETTINGS (LPM)
0-25 LPM	0-.5-1-2-3-4-5-6-8-10-15-25

## Specifications:

### Material:

Body Chrome-plated Brass  
 Internal Part: Brass core  
 T-handle: Chrome Plated Brass  
 Hose barb/DISS: Chrome Plated Brass

**Length:** Standard 5-5/8" (14.3 cm)  
**Weight:** Standard 16.4 oz. (465 grams)  
**Safety Relief Pressure:** 75-130 PSI  
**Diameter:** 1-3/8" (3.5 cm)  
**Maximum inlet pressure:** 3,000 psig (20.7 MPa)  
**Outlet pressure:** 50 psig (345 kPa)  
**Warranty:** 5 Years (see warranty statement for details)



**Genstar Technologies Co., Inc.**  
 4525 Edison Avenue  
 Chino, CA 91710 USA  
 TEL: 909.606.2726  
 FAX: 909.606.6485  
[www.gentehealthcare.com](http://www.gentehealthcare.com)